

Bruker Fourier300 NMR Spectrometer

Sample solution height must be 4.5 cm with sample fully dissolved in a deuterated solvent.

Open Bruker TopSpin 3.2 software. (Icon on desktop.)

Defining Sample


- Type **“new”** and enter information for your sample:
 NAME – **Folder name** (Your last name would be good). This will create a new folder for your dataset, or you may add to an existing folder by typing its name here.
 EXPNO – Increment this number for each sample you run.
 PROCNO – This value must remain “1”.
 Experiment – for ^1H , use **“PROTON”**; for ^{13}C , use **“C13CPD”**; for other options click “Select”.
 Set solvent – choose appropriate solvent from drop down menu.
 Leave “Execute getprsol” and “Keep Parameters” unchecked. Leave “Receivers” set to 1.
 DIR- Make sure the directory is set to “C:\Bruker\TopSpin3.2\data\nmrsu\nmr\Chem250”
 TITLE- Enter experimental details or the **sample name**.
- Click **OK**.

Inserting Sample (Check that dust cap has been removed.)

- Type **“ej”** to eject previous sample from NMR. Remove old sample from spinner, and carefully place your sample into the spinner. Use plastic depth gauge to check appropriate solution height: push tube through spinner gently to the bottom of gauge. Put your spinner and tube into the NMR sample probe. It should float.
- Type **“ij”** to lower the sample into the probe. You will hear a click when sample is properly positioned in probe.

Locking and Shimming


- Type **“lock”**. The “Solvents Table” window will pop up. Choose your solvent to lock on.
Watch progress at bottom left. Wait for “lockn: finished” but message disappears if move cursor.



- Type **“ro on”** to spin the sample. *Watch progress at bottom left.*
- Option: type **“gradshim”** or click on “Start Gradient Shimming”. *Watch progress at bottom.*
- Type **“getprosol”** to get standard acquisition parameters. *Watch progress at bottom left.*
- Option: Type **“ased”** to modify parameters such as number of scans (ns) or delay (D1).
- Type **“rga”** to set the receiver gain.

Collecting Data

- Type **“xaua”**. *Watch progress at bottom left.*
- Click on the ‘Acquire’ tab in the TopSpin menu bar to see the FID.



- You can check if you have enough data by typing **“tr”** and then processing as listed below. If enough samples have been acquired, type **‘halt’** to save all scans acquired and stop acquisition.

Processing Spectra

- Click on the ‘Process’ tab in the TopSpin menu bar.
- Select ‘Proc. Spectrum’ button to load, Fourier transform, and phase your data or type **“efp”**, then **“apk”**, then **“abs”**. For manual phase adjustment so that peaks are symmetric and go up, type **“ph”**. Right click at one end of the spectrum and select ‘Set Pivot Point’. Click and drag the ‘0’ in the upper left to adjust peaks near the pivot point; click and drag the ‘1’ to adjust the rest.
- Type **“pp”** or select ‘Pick Peaks’ button. Drag a box that includes the top of the peaks. Click on the Save+Return icon to keep values or return icon to exit. Type **“sref”** to set a specific value.
- Type **“.int”** or select ‘Integrate’ button. With the field goal icon yellow, click and drag over a peak to integrate. To set a specific integral value, right click and select ‘Calibrate Current Integral’.
- Click ‘Plot’ on the spectrum toolbar or type **“plot”**. Set paper settings to ‘Letter’. A useful layout file is “Standard KLB.xwp”. Zoom range can be specified by typing **“.zx”** or **“.all”**.

Always store a sample in spectrometer when finished! Type **“ro off”** and replace dustcap.

